



Powerology

**Full-Body
Massage Gun**

SKU: PSM002

Table of Contents

Overview	2
Adjusting the Angle	2
Product Description	3
Application	4
Specifications	4
Disclaimer	5
Disposal	5
Warranty	6
Contact Us	6

Overview

Created with you in mind, the 180° delivers results that everyone will benefit from. Whether your goal is relief from workout fatigue or chronic joint pain, to improve athletic performance and blood circulation, or simply to enjoy a relaxing massage, percussive massage therapy is the perfect solution.

Adjusting the Angle



Note:

1. Turn off the massager before adjusting the arm angle.
2. Press and hold the arm adjustment button, adjust the massage arm to the desired position, and then release the button.
3. There are three angle adjustment options: 45°, 0°, and 90°.

Product Description

Multilevel percussive frequencies penetrate deeply into your muscles and fascia (connective tissue) to provide long-lasting pain relief, decrease tension, improve blood circulation, and reduce inflammation. Equipped with six head attachments that target specific areas of your body, the 180° offers the precision you need.

- 1. Elevate your workouts:** Reduce muscle tension, improve flexibility, and decrease your chances of injury.
- 2. Recover faster and better:** Decrease inflammation, loosen tight joints, and break down scar tissue.
- 3. Relieve your pain:** Effectively treat trigger points, sore muscles, and chronic pain caused by tendinitis, bursitis, carpal tunnel, and more.
- 4. Relax more:** Boost your circulation and immune function, reduce stress levels, and get better sleep.

Note: Please consult with your doctor or physician before using this device if you have any injuries or medical conditions that may be of concern. This is not a medical device.



Application

The massager is a true game-changer in rehabilitation and sports performance. Massage guns have been shown to significantly enhance athletic performance and alleviate muscle soreness. The benefits are universal—from office and labor workers to professional athletes, weightlifters, and endurance athletes. Take full advantage of the versatility this product offers to achieve complete body recovery, helping you prepare better for your next workout.

Specifications

Battery Capacity	2000mAh
Input Voltage	5V
Working Time	3.5H
Charging Time	4.5H
Charging Port	Type-C
Speed Range	1500 to 3000RPM
Net Weight	920g
Product Size	22.5x17x10cm
Model Number	SM02

Disclaimer

- 1.** Read the instructions carefully before using the massager. The use of this equipment is subject to closely following the provided instructions. Our company is not responsible for any damage or loss resulting from inappropriate use of the device.
- 2.** Do not use the massager without the approval of your doctor if you are: pregnant, diabetic, have a pacemaker, recovering from surgery, suffer from migraines, neuropathy disorders, disc protrusion, or spinal inflammation.
- 3.** If you experience pain or discomfort while using the massager, stop immediately and consult your doctor or physiotherapist.

Disposal

This product must not be disposed of as unsorted household waste. It is important to separate such waste for proper treatment and recycling, in compliance with local waste management regulations.



Warranty

Products that you buy directly from our **Powerology** website or shop come with a 24-month warranty.

The 24-month warranty applies to products purchased directly from our **Powerology** website or store. If **Powerology** products are bought from any of our verified retailers, then the product is eligible for only a 12-month warranty. To extend your product's warranty, visit our website powerology.me/warranty and fill in your details in the provided form along with an uploaded picture of the product to process your request. Once approved, you will receive a confirmation email of the extended product warranty. Upload the required information within 48 hours of purchase to be eligible for a 24-month warranty period.

For more info, please check:
powerology.me/warranty

Contact Us

If you have any questions about this Privacy Policy, please contact us at: hey@powerology.me

Website: powerology.me

Instagram: [powerology_official](https://www.instagram.com/powerology_official)

Facebook: [powerology.ME](https://www.facebook.com/powerology.ME)